

# AVGPREP - Program for Computer Configuration Cloning Support

Document Revision 75.2 (8.2.2007)

**Copyright GRISOFT, s.r.o. All rights reserved.**

All other trademarks are property of their respective owners.

## Contents

<b>1. AVGPREP Program .....</b>	<b>3</b>
<b>2. The Process of Cloning.....</b>	<b>4</b>
2.1. Setting the default AVG Anti-Virus configuration .....	4
2.2. Using the AVGPREP program.....	4
2.3. Using the Sysprep Tool .....	4

## 1. AVGPREP Program

**AVGPREP** is a program intended primarily for administrators of large networks who may need to clone the computer's preset configuration to many PCs at once. This can be performed with use of the *Sysprep* tool.

The **AVG** license number is restricted to one computer only. During cloning the saved license data saved is devalued, and the **AVG** on the target stations reports the license number to be invalid.

To avoid this problem use the **AVGPREP** program. Before you start cloning the default configuration of the source computer with the help of *Sysprep*, it is necessary to use the **AVGPREP** program that creates a temporary *avgprep.dat* file, and with its use secures the automatic license correction on the target stations.

## 2. The Process of Cloning

If you want to clone the configuration of a computer to more PCs without re-registering every **AVG** installation on the target stations, follow these steps:

### 2.1. Setting the default AVG Anti-Virus configuration

Set the complete configuration intended for cloning on the source computer. The target stations will be installed with exactly the same settings, and the **AVG** license number used for the source station installation will be copied to all target stations!

### 2.2. Using the AVGPREP program

The **AVGPREP** program is available on the **AVG** installation CD, or you can download it from the Grisoft website at <http://www.grisoft.com>, Download section.

To launch **AVGPREP** use the program's executable file. You can launch it either directly with the specific parameter (the last four characters of the license number) or in the interactive mode where you will be invited to specify the last four characters of the license number.

Run the **AVGPREP** program to create the *avgprep.dat* temporary file into the **Avg7Data** directory. The **AVGPREP** program functionality is conditioned by the **AVG** installation (with the valid license number) on the source computer.

To avoid the **AVGPREP** program misapplying, it is necessary to enter the last four characters of the current license number (or the complete license number) to be able to create the *avgprep.dat* file. The generated file is then used in fully automatic mode, and the administrator does not have to deal with it any more.

### 2.3. Using the Sysprep Tool

With the help of **Sysprep** (or an alternative tool) create the source computer image, and clone its configuration to the target stations. Once the image is applied to the target computers, the **ProductId** changes, and the **AVG** licenses will be automatically corrected using the *avgprep.dat* file, and the file itself will be removed.